



INSTRUMENT SLEEVES


Instrument protection is a vital step in the reprocessing cycle. Using Key Surgical instrument protection, such as Instrument Sleeves, supports the protection of your valuable instruments from damage, as well as your team during handling. Use with hinged or non-hinged instruments to protect instrumentation during sterilization to point of use in the operating room.

“Ratcheted instruments should be unlatched. Racks, pins, stringers, or other specifically designed devices can be used to hold the instruments in the unlatched position.”

- AAMI/ANSI ST79, Section 8.2

ABOUT:

- Helps protect distal ends of instruments (hinged or non-hinged) from damage during sterilization and transportation
- Helps protect personnel from sharps injuries during handling
- Transparent pocket makes instrument visualization quick, safe, and clear
- Hinged instruments need to be in an open position to allow for proper sterilant contact of box locks and hinges during sterilization. Use instrument sleeve (01457) to keep hinged instruments open.
- Material is compatible with the following sterilization methods:
 - Steam
- Single use



PRODUCT NO.	DESCRIPTION	QUANTITY
01455	Instrument Protectors, paper/plastic, 2" x 5", fits instruments up to 32mm wide Designed for non-hinged instruments	100/pkg
01457	Instrument Protectors, paper/plastic, 3.5" x 6.6", fits instruments up to 40mm wide Designed for hinged instruments	100/pkg



a STERIS company

FOR MORE INFORMATION CALL 800.541.7995 OR VISIT [KEYSURGICAL.COM](https://www.keysurgical.com)

