




INSTRUMENT STRINGERS



Key Surgical® Instrument Stringers are ideal for organizing ring-handled instruments within a tray. Various widths are available to ensure that hinged instruments remain in an open position for effective washing and sterilization. Choose from several styles to fit your needs in the department.

- Durable, corrosion-resistant stainless steel (304 18/8 grade)
- Keeps ring-handled/hinged instruments in an open position for effective washing and sterilization
- Helps keep instruments organized within a tray
- Available widths: 2.5", 2.75", 3", 4", or 5.5"
- Available lengths: 4", 6", 8", 10", 12", 14", 16", 18", 20"
- Styles:
 - Ball & Socket
 - Pull-to-Open
 - Center Lock
 - Spring Lock
 - U-Shaped
 - T-Bar
- Compatible with steam sterilization. See IFU



"Ratcheted instruments should be unlatched. Racks, pins, stringers, or other specifically designed devices can be used to hold the instruments in the unlatched position."

-ANSI/AAMI ST79:2017, Section 8.2



**CONTACT YOUR STERIS REPRESENTATIVE
FOR MORE INFORMATION**

steris.com
Mentor, OH 44060 | USA

Key Surgical is a registered trademark of STERIS
MKTG4016 Rev G, 2023
M12628EN, 2023-11, Rev. A
©2023 STERIS. All rights reserved.

