



## TRACKING... DOWN TO THE INSTRUMENT LEVEL

Tracking inventory down to the instrument level doesn't have to be complicated. KeyDot is an effective, affordable asset tracking solution. This 2D data matrix barcode label is easy to apply, doesn't require instruments to be sent off-site, and with sterilization validation you can have peace of mind that choosing KeyDot will be a great choice for your investment.




- 2D data matrix barcode label designed to be used with compatible instrument tracking software\*

- Available in 4 sizes. (actual size shown left)
- Apply on stainless steel instruments or engineered thermo-set polymer (medical-grade plastic)

- Compatible with steam, VHP, and gas plasma sterilization and high-level disinfectants (HLDs)




1  
clean application area



2  
lift KeyDot from sheet



3  
place KeyDot on instrument



4  
apply firm pressure and sterilize to cure adhesive

## \*IS KEYDOT COMPATIBLE WITH YOUR SOFTWARE?

KeyDot is compatible with many software systems used throughout the industry. We suggest SPM® [Surgical Asset Tracking Software](#) from STERIS.

Contact your sales representative to learn more.

VISIT [KEYSURGICAL.COM](https://www.keysurgical.com) OR CONTACT YOUR SALES REPRESENTATIVE FOR MORE INFORMATION