

# CHANNEL CLEANING BRUSHES


Lumen Guard™



Surgical instrument cleaning just got a little smoother. Key Surgical® Lumen Guard Cleaning Brushes feature a smooth, inseparable, protective tip. Traditional exposed wire at the end of the brush is transformed through the use of Z-Tip™ technology, helping protect against potential scratches or damage to the lumen/channel.

- Channel cleaning brushes featuring **rigid nylon**, antimicrobial bristles on twisted stainless steel shafts with protective tip and loop end
- **Smooth, inseparable tip** protects against potential scratches inside lumen/channel
- Brush diameters range from .05" (1.27mm) to 1" (25.4mm)
- Lengths range from 6" to 24"
- Choose from 2/pkg or 50/pkg to **optimize inventory**

[See reverse side for individual product details](#)



VISIT [STERIS.COM/HEALTHCARE](https://www.steris.com/healthcare) OR CONTACT YOUR STERIS REPRESENTATIVE DIRECTLY

[steris.com](https://www.steris.com)  
Mentor, OH 44060 | USA

Key Surgical® is a registered trademark of STERIS  
MKTG4039 Rev D, 2024  
M12673EN. 2024-06, Rev. A  
©2024 STERIS. All rights reserved.

 **STERIS**

Product No.	Diameter (IN)	Diameter (MM)	FR	Bristle Area Length	Overall Length	Quantity
BZ10050 BZ1005050	.050"	1.27mm	4	1"	10"	2 50
BZ12062 BZ1206250	.062"	1.58mm	5	1"	12"	2 50
BZ18062 BZ1806250	.062"	1.58mm	5	2"	18"	2 50
BZ18079 BZ1807950	.079"	2.00mm	6	2"	18"	2 50
BZ12087 BZ1208750	.087"	2.20mm	7	1"	12"	2 50
BZ12098 BZ1209850	.098"	2.50mm	8	1"	12"	2 50
BZ18105 BZ1810550	.105"	2.67mm	8	2"	18"	2 50
BZ08125 BZ0812550	.125"	3.18mm	10	.75"	8"	2 50
BZ12125 BZ1212550	.125"	3.18mm	10	1"	12"	2 50
BZ18125 BZ1812550	.125"	3.18mm	10	2"	18"	2 50
BZ24125 BZ2412550	.125"	3.18mm	10	2"	24"	2 50
BZ18131 BZ1813150	.131"	3.33mm	10	2"	18"	2 50
BZ18144 BZ1814450	.144"	3.37mm	11	2"	18"	2 50
BZ12158 BZ1215850	.158"	4.00mm	12	1"	12"	2 50
BZ18158 BZ1815850	.158"	4.00mm	12	2"	18"	2 50
BZ18197 BZ1819750	.197"	5.00mm	15	2"	18"	2 50
BZ18236 BZ1823650	.236"	6.00mm	18	2"	18"	2 50
BZ16250 BZ1625050	.250"	6.35mm	19	2"	16"	2 50
BZ18250 BZ1825050	.250"	6.35mm	19	2"	18"	2 50
BZ18312 BZ1831250	.312"	8.00mm	24	2"	18"	2 50
BZ18354 BZ1835450	.354"	9.00mm	28	2"	18"	2 50
BZ18393 BZ1839350	.393"	10.00mm	30	2"	18"	2 50
BZ18433 BZ1843350	.433"	11.00mm	34	2"	18"	2 50
BZ18472 BZ1847250	.472"	12.00mm	N/A	2"	18"	2 50
BZ181000 BZ18100050	1"	25.40mm	N/A	2"	18"	2 50
BZ181543 BZ18154350	.157" to .433"	4mm to 11mm	N/A	2"	18"	2 50
BZ062550 BZ06255050	.250" to .500"	6.35mm to 12.7mm	N/A	.75"	6"	2 50