



POLYPECTOMY SNARES

Polypectomy Snares are used endoscopically in the removal of diminutive polyps, sessile polyps, pedunculated polyps and tissue from within the GI tract. Available in multiple styles (oval, mini-oval, or hexagonal) to suit your procedural needs.


ABOUT:

- Self-contained metal loop that is opened over the polyp and then closed to capture and remove the polyp tissue
- Available in the most preferred styles: oval and hexagonal
- Braided wire combines the strength of a firm snare with the flexibility of a soft snare to help with accurately capturing tissue. This combination can also help reduce the need for additional snare inventory in the GI suite
- Sterile
- Single Use



PRODUCT NO.	DESCRIPTION	WIDTH	QUANTITY
EN10300B	Stiff Multi-Wire, Hexagonal	10mm	10/box
EN10301B	Stiff Multi-Wire, Hexagonal	15mm	10/box
EN10303B	Stiff Multi-Wire, Hexagonal	25mm	10/box
EN10306B	Stiff Multi-Wire, Oval	10mm	10/box
EN10307B	Stiff Multi-Wire, Oval	15mm	10/box
EN10309B	Stiff Multi-Wire, Oval	25mm	10/box
EN10310B	Stiff Multi-Wire, Oval	30mm	10/box
EN 10330	Cold-Hot, Stiff Multi-Wire, Mini-Oval	10mm	10/box

YOU MIGHT ALSO LIKE



VISIT OUR WEBSITE KEYSURGICAL.COM OR CALL 800.541.7995 FOR MORE INFORMATION

