



POST-MORTEM BAGS

Post-mortem bags are intended for containment of a deceased human body. Featuring a c-bow opening for easy loading and high-strength reinforced material for peace of mind against tearing or ripping. One size fits most.

ABOUT:

- 35.4" W x 88.5" L (seam to seam)
- C-bow shaped zipper opening
- Intended for containment use only, not designed for transportation
- Heat sealed edges allow for leak resistance and strength

MATERIALS:

- Film: tear-resistant, reinforced LDPE (white) 0.2mm thickness
- Zipper Track: Nylon and Polyester
- Zipper: Zinc Alloy
- Contains no PVC
- Not made with natural rubber latex

PRODUCT NO.	QUANTITY
-------------	----------

M99925	10/box
--------	--------

