

# ROBOTIC ARM INSTRUMENT CARD



- The Robotic Arm Instrument Card is a **positioning tool** that holds the robotic arm instrument in a sterile packaging product (peel pouch or wrap)
- Designed for use **in conjunction with a sterile barrier** such as a peel pouch or sterilization wrap, when preparing for sterilization
- **Patented design** includes a viewing window for clear **visibility** of the distal tip
- Dedicated **slot** for placing **steam indicator**
- Features retrieval rings (one on each end) for aseptic presentation in the OR
- Compatible with steam sterilization (please refer to IFU)
- Single use



Product No.	Description	Qty
CSRC100	Robotic Arm Instrument Card	25

## RELATED PRODUCTS:

Product No.	Description	Qty
881091	Vis-U-All High Temp Tubing, 9"x 100'	6
PCC064	VERIFY Steam Integrating Indicator 4"	250
PCC065	VERIFY Steam Integrating Indicator 4"	1000
PCC067	VERIFY Steam Integrating Indicator 2"	500

**CONTACT YOUR STERIS REPRESENTATIVE TO LEARN MORE**