

SAFE CORD™

FLOOR STRIP



- Safe Cord is a durable, single-use safety strip designed to cover cables and cords on operating room floors to reduce the risk of staff injuries due to tripping hazards
- Tear-resistant material for durability during use and water resistant to prevent lifting of adhesive strips on the back
- Bright yellow color for easy visualization on the floor
- Support efficient O.R. turnover with ease of use and disposal
- Single-use. Not made with natural rubber latex



Product No.	Description	Qty
SC100	Safe Cord, 28" x 9", adhesive strips on back	30/box

**CONTACT YOUR STERIS REPRESENTATIVE
TO LEARN MORE**

steris.com
Mentor, OH 44060 | USA

MKTG4140 Rev A, 2023
Safe Cord™ is a registered trademark of Flagship Surgical
M13182EN.2023-12, Rev.A
©2023 STERIS. All rights reserved.

