



# TAGS: SCOPE REPROCESSING

Key Surgical® Scope Tags are designed to be used during scope reprocessing to record information such as scope model, flush, leak test pass/fail, and processing dates. The tags are double-sided, self-looping, and tear-resistant making them a perfect match for your communication needs in the department. Two different styles available to meet your various needs in scope reprocessing.


## SCOPE TAGS:

- Available in seven bright colors for easy color-coding
- Month, day, and date selection on one side
- 13.25" x 1.5"
- Self-looping tag
- Compatible with steam and gas plasma sterilization
- 100 tags/pkg
- Includes 1 lab marker



## SCOPE PROCESSING TAG:

- Capture detailed information (leak test, scope model, etc).
- 10.5" x 3"
- Self-looping tag
- Dedicated barcode placement area
- Compatible with steam and gas plasma sterilization
- 100 tags/pkg
- Includes 1 lab marker



a STERIS company