



# SHROUDGUARD™ COLUMN PROTECTORS

Prevent costly damage to the columns of STERIS® O.R. surgical tables by using ShroudGuard Column Protectors. Surgical table columns are often damaged due to inappropriately placed OR equipment (arm boards, sequential compression devices, convective warming devices, etc). The ShroudGuard Protector securely mounts to the surgical table, creating a secure fence around the column. It has been said that, “An ounce of prevention is worth... not paying thousands of dollars to repair the columns on your surgical table.”

## ABOUT:

- Rests on the table base cover, no modification to table required
- Designed for single person installation
- Easy to install and remove for cleaning and maintenance
- Keep your surgical tables in service and maximize OR uptime

PRODUCT NO.	DESCRIPTION	QTY
BF067	Column Protector for 3000-series STERIS/AMSCO Tables	1/box
BF068	Column Protector for 4000-series STERIS Tables	1/box
BF069	Column Protector for 5000-series STERIS Tables	1/box

## MATERIALS:

- Durable stainless steel construction
- Rubber bumpers prevent against scratches



BF067



BF068



BF069



a STERIS company

VISIT OUR WEBSITE [KEYSURGICAL.COM](http://KEYSURGICAL.COM) OR CALL 800.541.7995 FOR MORE INFORMATION

