




ENDOSCOPE MANUAL CLEANING: SMALL WONDER® CLEANING BRUSH

The Small Wonder cleaning brush is a specialized brush for the smallest and most delicate endoscope channels. Suited for cleaning the balloon inflation channel of EUS scopes, tiny channels of select bronchoscopes, choledochoscopes, and other small instrument channels ranging in size from 1.0mm to 2.25mm.

FOR CHANNEL
SIZES
1.0MM – 2.25MM



ABOUT:

- Available in lengths of 160cm or 230cm
- Flexible nylon bristles set on a non-slip, braided wire shaft to help with a steady grip
- Stainless steel cap at the end of the bristles is designed to help protect against scratches or damage as the brush moves through the channel
- Barbell handle at the end of the brush with grooves for grip

PRODUCT NO.	DESCRIPTION	BRISTLE DIAMETER (MM)	BRISTLE AREA LENGTH (MM)	SHEATH DIAMETER (MM)	CHANNEL DIAMETER (MM)	OVERALL LENGTH (CM)	QUANTITY
BX00711616	Small Wonder Cleaning Brush	2.25	6.25	0.5	1.0 - 2.25	160	15/box
BX00711617	Small Wonder Cleaning Brush	2.25	6.25	0.5	1.0 - 2.25	230	15/box



a STERIS company

VISIT OUR WEBSITE KEYSURGICAL.COM OR CALL 800.541.7995 FOR MORE INFORMATION

