



SOAKER SHEETS

Key Surgical® Soaker Sheets are specifically designed for heavy duty use, such as reprocessing orthopedic sets which can be the most challenging for packaging integrity and wet packs due to weight).

“Absorbent tray liners can absorb condensation and reduce the incidents of wet packs.”

-ANSI/AAMI ST79:2017

ABOUT:

- Use under trays to protect wraps from sharp corners and uneven surfaces of trays and baskets
- Designed for use with heavier trays (such as orthopedic trays) that are more prone to condensation challenges
- Provides rapid absorption of condensate during sterilization to help prevent wet packs
- May also be used as a shelf liner for drying racks
- Available in various sizes
- Non-linting, durable material
- Manufactured with low chloride and sulfate content
- Free from optical dyes
- Compatible with steam and EtO sterilization
- Single use



PRODUCT NO.	DESCRIPTION	DIMENSIONS	QUANTITY
SOK006A	Soaker Sheet, cellulose, 140gsm	250mm x 250mm (9.84" x 9.84")	800
SOK003A	Soaker Sheet, cellulose, 140gsm	305mm x 610mm (12" x 24")	400
SOK002A	Soaker Sheet, cellulose, 140gsm	500mm x 375mm (19.70" x 14.76")	200
SOK010A	Soaker Sheet, cellulose, 140gsm	500mm x 914mm (19.70" x 36")	96
SOK008A	Soaker Sheet, cellulose, 140gsm	560mm x 1219mm (22" x 48")	100
SOK009A	Soaker Sheet, cellulose, 140gsm	560mm x 1524mm (22" x 60")	100
SOK001A	Soaker Sheet, cellulose, 140gsm	750mm x 500mm (29.53" x 19.70")	96

VISIT [KEYSURGICAL.COM](https://www.keysurgical.com) OR CALL 800.541.7995 FOR MORE INFORMATION