




STERILIZATION WRAP

Strength. Flexibility. Peace of mind. Key Surgical sterilization wrap supports patient safety where it starts: sterilization. Our barrier wrap is strong (double-bonded SMS) to help protect against microorganisms (bacteria, viruses, spores) while allowing sterilization. The material, resistant to tears/rips, is also flexible so it's easy to work with and the unique pattern makes visualization for holes during aseptic presentation in the OR a breeze.

ABOUT:

- Up to **39% stronger*** material is paired with flexibility for a top-performing sterile barrier
- Each sheet is **tear-resistant**, double bonded (SMS)
- **Easy to use**, flexible wrap for securely wrapping trays, devices, and other surgical items
- Unique **diamond pattern** and **dual-color (blue/green)** supports visualization of holes or tears during aseptic presentation in the operating room
- Validated for steam, EtO, gas plasma, and vaporized hydrogen peroxide sterilization efficacy
- Single use

FOR WRAPPING TECHNIQUES, PLEASE REFER TO AAMI ST79



100 WEIGHT WRAP

PRODUCT NO.	DESCRIPTION	QUANTITY
WR100-24-24	24" x 24"	240/case

200 WEIGHT WRAP

PRODUCT NO.	DESCRIPTION	QUANTITY
WR200-18-18	18" x 18"	480/case
WR200-24-24	24" x 24"	240/case
WR200-30-30	30" x 30"	144/case
WR200-36-36	36" x 36"	144/case

300 WEIGHT WRAP

PRODUCT NO.	DESCRIPTION	QUANTITY
WR300-24-24	24" x 24"	240/case
WR300-30-30	30" x 30"	144/case
WR300-36-36	36" x 36"	144/case
WR300-54-54	54" x 54"	48/case

400 WEIGHT WRAP

PRODUCT NO.	DESCRIPTION	QUANTITY
S40024	24" x 24"	120/case
S40036	36" x 36"	72/case
WR400-40-40	40" x 40"	72/case
WR400-45-45	45" x 45"	48/case
S40048	48" x 48"	48/case

500 WEIGHT WRAP

PRODUCT NO.	DESCRIPTION	QUANTITY
WR500-18-18	18" x 18"	288/case
WR500-24-24	24" x 24"	120/case
WR500-30-30	30" x 30"	120/case
WR500-36-36	36" x 36"	72/case
WR500-45-45	45" x 45"	48/case
WR500-48-48	48" x 48"	48/case
WR500-54-54	54" x 54"	48/case

600 WEIGHT WRAP

PRODUCT NO.	DESCRIPTION	QUANTITY
S60036	36" x 36"	72/case
WR600-40-40	40" x 40"	48/case
WR600-45-45	45" x 45"	48/case
S60048	48" x 48"	48/case
S60054	54" x 54"	48/case



a STERIS company

VISIT OUR WEBSITE KEYSURGICAL.COM OR CALL 800.541.7995 FOR MORE INFORMATION



*ISO 9073-4 test method compared 100 weight wrap with market leading product showed a 39% improvement in tensile strength