



SURGICAL LIGHT HANDLE COVERS

STERIS surgical light handle covers, including disposable surgical light handle covers and surgical camera covers, are designed to fit over surgical light handles to provide sterile surgeon control and eliminate the need for cleaning and sterilizing of handles for every case. Designed for use with STERIS surgical lights and medical examination lights. All styles are sterile, disposable, and not made with natural rubber latex.

ABOUT:

Compatible with select STERIS Harmony and HarmonyAIR Surgical Lights including the Harmony LA300 Exam Light. The covers provide sterile control of power on/off, intensity levels and pattern size. Individually packaged.

PRODUCT NO. QUANTITY

LB53 Single pack; (120 unit count)



ABOUT:

Compatible with select STERIS Harmony and HarmonyAIR Surgical Lights including the Harmony LA300 Exam Light with high definition in-light cameras. The covers provide sterile control of power on/off, intensity levels and pattern size. Individually packaged.

PRODUCT NO. QUANTITY

LB54 SGL pack; 30/per case



ABOUT:

Universal Light Handle Cover is compatible for use with universal surgical light handle adapters. Dual packaged.

PRODUCT NO. QUANTITY

LB82 DL pack 200/bx

