




SWINGING BIOPSY FORCEPS

Obtain biopsy tissue from the gastrointestinal tract in an efficient manner with Swinging Biopsy Forceps from Key Surgical. The pivot feature of Swinging Biopsy Forceps helps reduce, if not eliminate, the need for additional maneuvering of the endoscope within the GI tract.

ABOUT:

- Swinging function helps improve physician's ability to capture tissue needed for biopsy
- Available in both alligator and oval cup designs (with or without the positioning spike)
- Length markers to assist with insertion and withdrawal process
- Ergonomic handle
- Fenestrated (allows view of tissue through the forceps)
- Coated – to help with insertion
- Compatible with 2.8mm biopsy channels (max. diam of 2.4mm/working lengths of 180cm/230cm)
- Sterile
- Single use



Alligator Cup Without Spike

Alligator Cup With Spike

Oval Cup Without Spike

Oval Cup With Spike



VISIT OUR WEBSITE [KEYSURGICAL.COM](https://www.keysurgical.com) OR CALL 800.541.7995 FOR MORE INFORMATION





SWINGING BIOPSY FORCEPS

PRODUCT NO.	DESCRIPTION	DIAMETER	LENGTH	QUANTITY
EN10267AS	Alligator cup without spike	2.4mm	1800mm	10/box
EN10268AS	Alligator cup without spike	2.4mm	2300mm	10/box
EN10269AS	Alligator cup with spike	2.4mm	1800mm	10/box
EN10270AS	Alligator cup with spike	2.4mm	2300mm	10/box
EN10271AS	Oval cup without spike	2.4mm	1800mm	10/box
EN10272AS	Oval cup without spike	2.4mm	2300mm	10/box
EN10273AS	Oval cup with spike	2.4mm	1800mm	10/box
EN10274AS	Oval cup with spike	2.4mm	2300mm	10/box



VISIT OUR WEBSITE [KEYSURGICAL.COM](https://www.keysurgical.com) OR CALL 800.541.7995 FOR MORE INFORMATION

