



# VASCULAR BOOTIES


Handling fine and delicate sutures can be challenging. Key Surgical® Vascular Booties are soft instrument covers designed to help effectively and firmly grasp suture material to prevent slipping. Available in two sizes and 5 colors for color-coding needs.

## ABOUT:

- Radiopaque
- BPA free
- Sold sterile, ready-to-use
- Packaged in pads (2, 3, or 5 pairs per pad)
- Adhesive-backed pads for universal placement
- Sold 5 packages per box
- Colors available: white, blue, yellow, red, pink, or assorted
- Not made with natural rubber latex
- Single use



## AVAILABLE COLORS



## ALSO AVAILABLE IN PAIRS OF 2 OR 3 (YELLOW, .062" x .437" ONLY)

