



# PERSONAL PROTECTIVE EQUIPMENT: VR® SLEEVE GLOVE


Protect your arms and hands during manual cleaning with a VR Sleeve Glove. This PPE is a unique fusion of a durable nitrile glove and a protective sleeve and at 28" long, provides complete hand and arm protection. You don't have to compromise comfort for durability with the VR Sleeve Glove; it's durable to protect against a broad range of chemicals, yet flexible so it works *with* you, not the other way around.

## ABOUT:

- Complete, durable hand and arm protection through a unique fusion of a 15 mil, lined nitrile glove and a 6 mil protective sleeve
- Overall length of 28"
- Comfortably flexible with high tactile sensitivity
- Elastic cuff for adjustable positioning on upper arm
- Meets or exceeds OSHA Level D chemical protective apparel requirements for splash protection
- Resistant to a broad range of chemicals including hydrofluoric, hydrochloric, and sulfuric & nitric acids
- Pairs are individually packed for quick and easy use
- Made in the USA
- Not made with natural rubber latex
- Available in sizes S, M, L, XL (*see size chart on back for glove sizing*)



Lined for comfort



PRODUCT NO.	DESCRIPTION	QUANTITY
PS41380	VR Sleeve Glove, 80 Series, Small (7)	100/case
PS41480	VR Sleeve Glove, 80 Series, Medium (8)	100/case
PS41580	VR Sleeve Glove, 80 Series, Large (9)	100/case
PS41680	VR Sleeve Glove, 80 Series, Extra Large, (10)	100/case




a STERIS company

FOR MORE INFORMATION CALL 800.541.7995 OR VISIT [KEYSURGICAL.COM](http://KEYSURGICAL.COM)



# TO DETERMINE YOUR GLOVE SIZE:

1. Print the illustration below at 100% size.
2. Place your right hand on the diagram.
3. The blue line should be between your thumb and index finger.



a STERIS company

FOR MORE INFORMATION CALL 800.541.7995 OR VISIT [KEYSURGICAL.COM](https://www.keysurgical.com)

