



STERILIZATION RACKS

Sterilization racks designed to help you conveniently wash, sterilize, store and dispense various lengths and diameters of K-wires and pins. Made from anodized aluminum, the racks are lightweight but made to withstand repeated use. Dimensions clearly marked on the rack to assist with quick identification. Folds closed but easily converts to a table top stand for use in the O.R.

K-WIRE AND PIN RACK*

FOR USE WITH 9" K-WIRES AND PINS. HOLDS SIX OF EACH IN THESE DIAMETERS:

.028" (.7mm)	.054" (1.4mm)	.093" (2.4mm)	.140" (3.5mm)
.035" (.9mm)	.062" (1.6mm)	.109" (2.8mm)	.156" (4.0mm)
.045" (1.1mm)	.078" (2.0mm)	.125" (3.2mm)	.177" (4.5mm)



KI-71-RACK
12.5" x 5.25" x 1.5"

K-WIRE RACK*


FOR USE WITH 4" - 6" K-WIRES. HOLDS SIX OF EACH IN THESE DIAMETERS:

.028" (.7mm)	.045" (1.1mm)	.062" (1.6mm)
.035" (.9mm)	.054" (1.4mm)	



KI-71-RACK-46
6.25" x 4.125" x .75"

RELATED PRODUCTS:



K-wires and Steinmann Pins

Pin covers



FOR MORE INFORMATION CALL 800.541.7995 OR VISIT KEYSURGICAL.COM



*Pins and Wires not included