



STRINGER ASSISTANT

We can all agree that a little help goes a long way. This is especially true in the fast-paced environment of instrument reprocessing. The Key Surgical® Stringer Assistant is designed to help your team be efficient with instrument stringing. The Stringer Assistant can fit multiple stringer widths and with its anti-skid base, you will wonder how you ever got along without it.

“Ratcheted instruments should be unlatched. Rack, pins, stringers, or other specifically designed devices can be used to hold the instruments in the unlatched position.”

-ANSI/AAMI ST79:2017, Section 8.2


ABOUT:

- Solid acrylic base
- 10.25" x 7.5"
- Anti-skid bottom
- 7" slot for stringer (tapered)
- White
- 1/pkg



PRODUCT NO.	DESCRIPTION	DIMENSIONS
IS-675105	Acrylic	10.25" x 7.5"

YOU MIGHT ALSO LIKE:



Instrument Stringers

