




SURGICAL TRAY CORNER PROTECTORS

A rip in the sterilization wrap can certainly put a major wrinkle in your day. Protect wrapped sterilization trays from rips and tears with the use of surgical tray corner protectors and reduce unnecessary costs associated with contaminated trays. We offer disposable or reusable styles to meet your specific needs.

REUSABLE PROTECTORS:

- Medical grade silicone
- Lime green color (for better visibility against wrap)
- Compatible with steam and gas plasma sterilization
- Reusable (tested to withstand a minimum of 15 sterilization cycles)
- Not made with natural rubber latex
- 2.8" x 2.8" x 1.2"

PRODUCT NO. **QUANTITY**
 CS-9113 80/box



DISPOSABLE PROTECTORS:

- Foam material, blue
- Compatible with steam sterilization only
- Single use
- 6" x 6" x 1.25"

PRODUCT NO. **QUANTITY**
 CS-9125C 80 sheets/box
 (4 pre-cut squares/sheet)



YOU MIGHT ALSO LIKE:



KEYSURGICAL®