



# TRAY LINERS

Help guard against wet packs with Key Surgical® Tray Liners. Available in multiple sizes and with or without sides, these liners help facilitate the sterilization process by absorbing moisture within or under a tray.

“Absorbent tray liners can absorb condensation and reduce the incidents of wet packs.”

-ANSI/AAMI ST79:2017

## ABOUT:

- Helps absorb moisture during sterilization
- Available in multiple sizes
- Single use



TL-105-250



TL-101-200



TL-204-500



TL-202-500



TL-203-500

PRODUCT NO.	DESCRIPTION	QUANTITY
TL-105-250	7.09" (180mm) x 11.02" (280mm)	250/pkg
TL-101-200	11.81" (300mm) x 19.69" (500mm)	200/pkg
TL-204-500	13.78" (350mm) x 13.78" (350mm), with sides	500/pkg
TL-202-500	13.78" (350mm) x 22.83" (580mm), with sides	500/pkg
TL-203-500	15.75" (400mm) x 23.62" (600mm), with sides	500/pkg

